



# Non-Crystalline Diffraction and Scattering at the BioCAT Facility at the Advanced Photon Source




Tom Irving

BioCAT, CSRRI and Dept. BCPS, Illinois  
Institute of Technology, Chicago IL



# What is BioCAT?

- A NIH-supported research center for the study of partially ordered and disordered biological materials
- Comprises an undulator based beamline, (18-ID) associated laboratory and computational facilities.
- Available to all scientists on basis of peer-reviewed beamtime proposals





# Scientific Mission of BioCAT

---

Modality	Applications
Fiber Diffraction	Muscle, Connective Tissue, virus structure, DNA, Amyloids
Solution Scattering	Protein/RNA folding, Protein-ligand interactions
Scanning $\mu$ -x-ray florescence microscopy and diffraction imaging	Metal metabolism, neuro-degenerative disease, cancer

---

# SAXS Instrument on the BioCAT 18ID - Undulator Beamline



- Total X-ray flux  $1\text{-}2.5 \times 10^{13}$  photons/s
- Focal spot size ranges from  $< 50 \mu\text{m}$  vertical and  $< 150 \mu\text{m}$  horizontal to  $\sim 3 \times 1.5 \text{ mm}$



## **Why is Doubly Focused Undulator Radiation Good for SAXS/NCD Studies?**

- High flux density for time-resolved applications and obtaining good counting statistics from weakly scattering systems
- Small focal spots along with well-matched high resolution area detectors allow resolving weak diffraction features in the presence of high backgrounds typical of biological samples
- Spatially resolved systems (Micro-diffraction)

# Aviex PCCD 16080 Detector



**160 x 80 mm**

**2000 x 4000 39  $\mu$ m  
pixels**

**$\sim 5/\text{ADU's}/12 \text{ keV}$   
Photon**

**$\sim 1 \text{ ADU}$  read noise**

**$\sim 1 \text{ s}$  readout**

**Flexible binning modes**

**Both high sensitivity and high spatial resolution**

# Pilatus 100K Photon Counting Detector



2D detector, continuous readout  
till computer disk full

Pixel size             $172 \times 172 \mu\text{m}^2$   
Format    $487 \times 195 = 94\,965$  pixels  
Active area         $83.8 \times 33.5 \text{ mm}^2$   
Counting rate       $> 2 \times 10^6$  counts/s/pixel  
Energy range        $3 - 30 \text{ keV}$   
Readout time        $< 2.7 \text{ ms}$   
Framing rate        $> 200 \text{ Hz}$   
Power consumption  $5\text{W}$ , air cooled


Limited Energy  
resolving ability


200 ns  
electronic  
gating possible

**Stretch  
Activated**

**Relax**

**Low  
gain**






# Micro-Diffraction

Minimum focal spot size  $\sim 5 \times 5 \mu\text{m}$

Can adjust focal spot position on the fly to optimize flux vs. divergence

CCD detector with  $4000 \times 4000$   $40 \mu\text{m}$  pixels



Automated “Diffraction Mapping” possible






# Structural Studies of Collagen Type 2 from Lamprey Notochord(Orgel)

Histology mapping shows where the most crystalline /oriented regions are located



1.1 nm equatorial

# Solution Scattering



Camera lengths from

~ 1 m to 3.5 m

Standard lengths:



2 m  $q$  range:

$0.007 \text{ \AA}^{-1}$  to  $0.4 \text{ \AA}^{-1}$

3m  $q$  range:

$0.005 \text{ \AA}^{-1}$  to  $0.3 \text{ \AA}^{-1}$


# FPLC + SAXS



Plasminogen data courtesy N. Menhart, IIT




# Stopped Flow for Kinetics



- Bio-Logic SFM-400 stopped-flow with new Biologic MEC 22998 microvolume mixer
- $\sim 0.5$  ms dead time
- Requires flow rate of 8ml/s of sample.
- Typical sample consumption  $\sim 250 \mu\text{l}$  (mixed volume) per mixing event



# Stopped Flow with Pilatus Detector



- Continuous SAXS data acquisition with a time resolution of 1 ms with a reaction dead time of 0.5 ms.


Data courtesy of Briber/Woodson group

# Sub-Millisecond Time Resolved Studies with Turbulent Flow Mixers (Osman Bilsel)




4.5 M  $\rightarrow$  0.45 M GdnHCl refolding (with 0.2M imidazole)  
of cytochrome c

Time resolution :  $\sim 80$  microseconds



# BioCAT Staff

- **Dr. Olga Antipova**
- **Dr. Srinivas Chakravarthy**
- **Mr. Rich Heurich**
- **Prof. Tom Irving**
- **Ms. Clareen Krolik**
- **Prof. Joseph Orgel**
- **Mr. Mark Vukonich**
- **Dr. Weifeng Shang**



**BioCAT is supported by  
the : National Institute of  
General Medical Sciences  
9 P41 GM103622-17**

**Part of the US National  
Institutes of Health**